

HYGIENE

AUSTRIA



Mouth-Nose Mask, Type-I

- Made in Austria
- 2. Certified
- 3. 3 layers
- 4. High breathability due to high quality non-woven materials



1. Wash your hands using soap and water or a hand sanitizer.



2. Pull the loops of the mask back to fit them around your ears.



3. Lightly press the metal piece in your mask across your nose to fit.



4. Pull the upper and lower sides of your mask over nose and chin.



Avoid touching the mask while wearing.





Mask Specifications

| Brand | - Hygiene Austria LP GmbH |
|-----------------------------|---|
| Name | - Mouth-Nose Mask, Type-I |
| Туре | - Non-Medical |
| Quantity/Units | - 50 pc. per unit |
| Colour | - Blue-white, upon request |
| Instructions and Durability | - Shelf life: 5 years - When humid, the mask is to be removed and disposed of - Non-reusable, disposable |
| Storage Conditions | - Store above 0°C and below 40°C in a dry place without direct sunlight |
| Size | - 175 x 95 mm (slight deviations possible) |
| Weight | - 6,21 g (slight deviations possible) |
| Certificates | - EN 14683:2019, Type-I - CE certified machines - TÜV tested - OEKO-TEX® Standard 100, Certificate #VN020 168438 OETI |
| Product specifications | Made out of 3 layers of non-woven material, 2 outside layers spunbond nonwovens and 1 inner layer melt-blown nonwovens polypropylene 8cm rubber-coated, extra strong and flexible nose piece |

Outer Carton Specifications

| Measurements Carton | - 580 x 310 x 400 mm |
|------------------------|----------------------------|
| Quantity | - 36 units, 50 pieces each |
| Material | - Corrugated cardboard |





Packaging

Packaging available in 50 piece format.

Through a light perforation, the masks can be easily pulled out from the top of the packaging unit.



